








# Precept<sup>®</sup>



## CHILDREN'S FACE MASK

 **15150**

-  **CHILD-SIZED TO IMPROVE FIT AND PROTECTION FOR CHILDREN**
-  **COLORFUL FUN PRINT TO ENCOURAGE COMPLIANCE**
-  **EASY TO USE EARLOOPS AND INDIVIDUALLY FORMABLE NOSE WIRE FOR SECURE, COMFORTABLE FIT**
-  **SOFT, LIGHTWEIGHT FABRIC ENHANCES COMFORT**
-  **MEETS ASTM F2100 REQUIREMENTS FOR BFE, PFE,  $\Delta$ P AND FLAMMABILITY**
-  **LATEX FREE**



ITEM	DESCRIPTION	COLOR	PACKAGING
15150	Children's Face Mask	Happy Face Print	75/Bx, 10Bx/Cs

*As recommended by the Centers for Disease Control and Prevention (CDC), to decrease transmission of respiratory infections, Precept's Children's Face Mask is designed specifically to be used by Children ages 3 – 12, in healthcare facilities, when exhibiting symptoms of a respiratory infection.*