

PROCEDURE MASKS

Item #	Description	BFE (3 µm)	PFE (0.1 µm)	Delta P	Flammability
15100	Isolation Mask	≥98%	≥98%	Pass	Class 1
15101	Procedure Mask	≥98%	≥98%	Pass	Class 1
15110	Procedure Mask	≥99%	≥98%	Pass	Class 1
15111	Procedure Mask	≥99%	≥98%	Pass	Class 1
15120	Procedure Mask	≥98%	≥97%	Pass	Class 1
15130	KittyKare Procedure Mask	≥99%	≥99%	Pass	Class 1

SURGICAL MASKS

Item #	Description	BFE (3 µm)	PFE (0.1 µm)	Delta P	Flammability
15200	Lite & Cool Surgical Mask	≥98%	≥98%	Pass	Class 1
15201	Standard Surgical Mask - Classic Style	≥98%	≥98%	Pass	Class 1
15215	Sensitive Skin Surgical Mask	≥98%	≥98%	Pass	Class 1
15220	Lite Pouch Surgical Mask	≥99%	≥99%	Pass	Class 1
15230	KittyKare Surgical Mask	≥98%	≥98%	Pass	Class 1
15250	Lite & Cool Horizontal Tie Mask	≥99%	≥98%	Pass	Class 1
15600	Standard Surgical Mask	≥99%	≥98%	Pass	Class 1
65-3333	Comfort-Cone Surgical Mask	≥98%*	N/A	Pass	Class 1

*Modified Greene & Vesley Method

SPECIALTY MASKS

Item #	Description	BFE (3 µm)	PFE (0.1 µm)	Delta P	Flammability
15210	Foam Anti-Fog Surgical Mask	≥98%	≥98%	Pass	Class 1
15211	Film Anti-Fog Surgical Mask	≥98%	≥98%	Pass	Class 1
15212	Tape Anti-Fog Surgical Mask	≥98%	≥98%	Pass	Class 1
15216	Sensitive Skin Foam Anti-Fog Surgical Mask	≥98%	≥98%	Pass	Class 1

FLUID RESISTANT MASKS

Item #	Description	ASTM F 1862 (Blood Splash Resistance)	BFE (3 µm)	PFE (0.1 µm)	Delta P	Flammability
15300	FluidGard™ Procedure Mask	160 mm Hg	≥99%	≥99%	Pass	Class 1
15301	FluidGard™ Anti-Fog Procedure Mask	160 mm Hg	≥99%	≥99%	Pass	Class 1
15310	FluidGard™ Procedure Mask w/ Extended Shield	160 mm Hg	≥99%	≥99%	Pass	Class 1
15315	FluidGard™ Procedure Mask w/ SnapShield	160 mm Hg	≥99%	≥99%	Pass	Class 1
15320	FluidGard™ Surgical Mask	160 mm Hg	≥99%	≥99%	Pass	Class 1
15330	FluidGard™ Surgical Mask w/ Extended Shield	160 mm Hg	≥99%	≥99%	Pass	Class 1
15335	FluidGard™ Surgical Mask w/ SnapShield	160 mm Hg	≥99%	≥99%	Pass	Class 1
15340	FluidGard™ Procedure VeilMask	120 mm Hg	≥99%	≥99%	Pass	Class 1
15345	FluidGard™ Procedure VeilMask w/ Extended Shield	120 mm Hg	≥99%	≥99%	Pass	Class 1
15350	FluidGard™ Surgical VeilMask	120 mm Hg	≥99%	≥99%	Pass	Class 1
15355	FluidGard™ Surgical VeilMask w/ Extended Shield	120 mm Hg	≥99%	≥99%	Pass	Class 1
15400	FluidGard™ KittyKare Procedure Mask	160 mm Hg	≥99%	≥99%	Pass	Class 1
15410	FluidGard™ KittyKare Procedure Mask w/ Ext. Shield	160 mm Hg	≥99%	≥99%	Pass	Class 1
15420	FluidGard™ KittyKare Surgical Mask	160 mm Hg	≥99%	≥99%	Pass	Class 1
15500	Fluid Resistant Procedure Mask	160 mm Hg	≥99%	≥99%	Pass	Class 1
15510	Fluid Resistant Procedure Mask w/ Extended Shield	160 mm Hg	≥99%	≥99%	Pass	Class 1
15520	Fluid Resistant Surgical Mask	160 mm Hg	≥99%	≥99%	Pass	Class 1
15530	Fluid Resistant Surgical Mask w/ Extended Shield	160 mm Hg	≥99%	≥99%	Pass	Class 1
65-3395	N95 Particulate Respirator (Regular Size)	160 mm Hg	≥99%	≥99%	Pass	Class 1
65-3395s	N95 Particulate Respirator (Small Size)	160 mm Hg	≥99%	≥99%	Pass	Class 1

Face Masks



PROCEDURE MASKS



Precept's procedure masks combine protection with comfort, to ensure ultimate satisfaction. Masks are constructed of quality materials, including tissue, polypropylene, and a wet-laid nonwoven made from cellulose. Enclosed nosepieces are malleable to ensure a snug, comfortable fit.

SURGICAL MASKS

Precept's surgical masks offer comprehensive protection, with bacterial filtration efficiencies exceeding 98%. The inner mask material has excellent fiber pindown, to help prevent facial irritation. Masks are cool and breathable, excellent for use during long procedures.



SPECIALTY MASKS



The specialty mask line incorporates various types of anti-fog masks. Precept offers several anti-fog masks made with soft foam, including a mask for sensitive skin. Tape anti-fog masks feature non-allergenic, pressure sensitive tape that provides a secure seal. Precept's film anti-fog mask incorporates a plastic film strip to help ensure optical clarity.

N95 PARTICULATE RESPIRATORS

Precept's N95 Particulate Respirators have been approved by NIOSH (the National Institute for Occupational Safety & Health). These N95 respirators are also highly fluid resistant, having passed the ASTM F 1862 blood splash test at the highest level of 160 mm Hg.



FLUID RESISTANT MASKS

Precept's fluid resistant masks provide the highest level of protection against splash from bodily fluids. These masks pass the ASTM F 1862 blood splash test at the highest level of 160 mm Hg. The masks offer excellent protection against bacteria, with bacterial filtration efficiencies exceeding 99%. Masks are available with or without anti-fog, non-glare extended shields. Fluid resistant masks come in a variety of prints, including blue, blue diamond, or Precept's exclusive KittyKare print. FluidGard masks are also available in Precept's innovative VeilMask design.



VEIL MASKS



VeilMasks offer excellent protection against bacteria and splash from bodily fluids. The fluid resistant veil material extends from the bottom of the face mask to the upper chest. This protects the neck area which is left exposed by traditional face masks. VeilMasks pass the ASTM F 1862 blood splash test at 120 mm Hg and filter against more than 99% of bacteria. VeilMasks are also available with and without extended shields.

FACE SHIELDS

Precept offers two sizes of face shields: 3/4 length and full length. The foam strip cushions against the forehead and the adjustable headband helps provide a secure, comfortable fit. Each shield is treated with an anti-fog, anti-glare coating to ensure optical clarity.

