



## Fluid Resistant Face Masks



15710



15320



15330



**Precept's Fluid Resistant Masks** provide the **highest level of protection** against splash from bodily fluids. These masks pass the ASTM 1862 Blood Splash Test at the **highest level of 160 mmHg**. They have also met ASTM 2100 Performance Requirements to be labeled as High Barrier. These masks offer **outstanding protection** against bacteria and particulates, with both bacterial filtration efficiencies (BFE) and particulate filtration efficiencies (PFE) exceeding 99%. Masks are available with or without anti-fog, non-glare extended shields, come in a variety of prints, and include **Precept's exclusive GlareGard™ Masks**, for the ultimate in glare reduction, and **Fluid Resistant, Anti-Fog Laser Mask**.

Item #	Description	Features	Color	Packaging
15300	FluidGuard® Procedure Mask	<ul style="list-style-type: none"> <li>■ Sonic welded for durability</li> </ul>	Blue Diamond	50/Bx, 10 Bx/Cs
15301	FluidGuard® Anti-Fog Procedure Mask	<ul style="list-style-type: none"> <li>■ Soft Foam Anti-Fog Strip increases comfort and inhibits eyewear fogging</li> </ul>	Blue Diamond	50/Bx, 10 Bx/Cs
15310	FluidGuard® Procedure Mask with Extended Shield	<ul style="list-style-type: none"> <li>■ Extended shield maximizes facial protection</li> <li>■ Soft Foam Anti-Fog Strip increases comfort and inhibits eyewear fogging</li> </ul>	Blue Diamond	25/Bx, 4 Bx/Cs
15330	FluidGuard® Surgical Mask with Extended Shield	<ul style="list-style-type: none"> <li>■ Black Anti-Glare Strip reduces reflections and glare</li> </ul>	Blue Diamond	25/Bx, 4 Bx/Cs
15315	FluidGuard® Procedure Mask with SnapShield	<ul style="list-style-type: none"> <li>■ Adjustable, removable Snap-on Shield maximizes facial protection</li> <li>■ Soft Foam Anti-Fog Strip increases comfort and inhibits eyewear fogging</li> <li>■ Black Anti-Glare Strip reduces reflections and glare</li> </ul>	Blue Diamond	25/Bx, 4 Bx/Cs
15320	FluidGuard® Surgical Mask	<ul style="list-style-type: none"> <li>■ Soft Foam Anti-Fog Strip increases comfort and inhibits eyewear fogging</li> <li>■ Sonic welded for durability</li> </ul>	Blue Diamond	50/Bx, 6 Bx/Cs
15335	FluidGuard® Surgical Mask with SnapShield	<ul style="list-style-type: none"> <li>■ Adjustable, removable Snap-on Shield maximizes facial protection</li> <li>■ Soft Foam Anti-Fog Strip increases comfort and inhibits eyewear fogging</li> <li>■ Black Anti-Glare Strip reduces reflections and glare</li> </ul>	Blue Diamond	25/Bx, 4 Bx/Cs
15400	FluidGuard® KittyKare™ Procedure Mask	<ul style="list-style-type: none"> <li>■ Child-friendly outer print of blue and pink kittens</li> <li>■ Sonic welded for durability</li> </ul>	KittyKare™ Print	50/Bx, 10 Bx/Cs

## Fluid Resistant Face Masks



Item #	Description	Features	Color	Packaging
15410	FluidGard <sup>®</sup> KittyKare <sup>™</sup> Procedure Mask with Extended Shield	<ul style="list-style-type: none"> <li>Extended shield maximizes facial protection</li> <li>Child-friendly outer print of blue and pink kittens</li> <li>Soft Foam Anti-Fog Strip increases comfort and inhibits eyewear fogging</li> <li>Black Anti-Glare Strip reduces reflections and glare</li> </ul>	KittyKare <sup>™</sup> Print	25/Bx, 4 Bx/Cs
15420	FluidGard <sup>®</sup> KittyKare <sup>™</sup> Surgical Mask	<ul style="list-style-type: none"> <li>Child-friendly blue and pink kitten print</li> <li>Sonic welded for durability</li> </ul>	KittyKare <sup>™</sup> Print	50/Bx, 6 Bx/Cs
15500	Fluid Resistant Procedure Mask	<ul style="list-style-type: none"> <li>Soft Foam Anti-Fog Strip increases comfort and inhibits eyewear fogging</li> </ul>	Blue	50/Bx, 10 Bx/Cs
15510	Fluid Resistant Procedure Mask with Extended Shield	<ul style="list-style-type: none"> <li>Extended shield maximizes facial protection</li> <li>Soft Foam Anti-Fog Strip increases comfort and inhibits eyewear fogging</li> </ul>	Blue	25/Bx, 4 Bx/Cs
15530	Fluid Resistant Surgical Mask with Extended Shield	<ul style="list-style-type: none"> <li>Black Anti-Glare Strip reduces reflections and glare</li> </ul>		
15520	Fluid Resistant Surgical Mask	<ul style="list-style-type: none"> <li>Soft Foam Anti-Fog Strip increases comfort and inhibits eyewear fogging</li> </ul>	Blue	50/Bx, 6 Bx/Cs
15525	Laser Mask, Fluid Resistant Anti-Fog	<ul style="list-style-type: none"> <li>Filters over 99% of particles as small as 0.1 micron</li> <li>Soft Foam Anti-Fog Strip Cushions Face &amp; Inhibits Fogging</li> <li>Contains No Dyes or Colorants on inner or outer facing or ties</li> </ul>	White	50/Bx, 6 Bx/Cs
15710	GlareGard <sup>™</sup> Procedure Mask with Extended Shield	<ul style="list-style-type: none"> <li>Extended Shield, excellent fluid resistance, and superior bacterial &amp; particulate filtration maximize protection</li> <li>Expanded Black Anti-Glare Strip, soft foam anti-fog strip, and anti-fog shield coating maximize visibility</li> </ul>	Black/Blue	25/Bx, 4 Bx/Cs
15730	GlareGard <sup>™</sup> Surgical Mask with Extended Shield	<ul style="list-style-type: none"> <li>Extra smooth inner fabric layer, ties, and bindings minimize facial irritation and ensure a comfortable fit</li> </ul>		
65-3345	DualGard <sup>™</sup> Procedure Mask with Extended Shield	<ul style="list-style-type: none"> <li>Adjustable, removable Snap-on Shield maximizes facial protection</li> <li>Anti-Fog Shield Coating enhances visibility</li> </ul>	Blue	25/Bx, 4 Bx/Cs
65-3347	DualGard <sup>™</sup> Procedure Mask with Foam Shield	<ul style="list-style-type: none"> <li>Adjustable, removable Snap-on Shield maximizes facial protection</li> <li>Soft Foam Anti-Fog Strip cushions forehead and inhibits eyewear fogging</li> </ul>	Blue	25/Bx, 4 Bx/Cs