

Now you see it!

Competitive Face Mask
with a Reflection



Now you don't!

GlareGard™ Face Mask
without a Reflection



Compare Precept's GlareGard™ Fluid Resistant Face Masks with your current Fluid Resistant Face Masks with Shields

Precept's GlareGard™ Masks with Extended Shield alleviate irritating reflections that adversely affect visibility, minimize glare, and provide excellent protection against bacteria and contamination from bodily fluids. GlareGard™ masks exceed 99% Bacterial and Particulate Filtration Efficiencies and the masks pass the ASTM F 1862 Blood Splash Resistance Test at the highest tested level (160 mm Hg). The soft foam strip and anti-fog coating help prevent the shield from fogging. Lightweight materials and soft bindings and ties ensure a snug, yet comfortable fit.

Item #	Description	Features	Color	Packaging
15730	GlareGard™ Surgical Mask with Extended Shield	<ul style="list-style-type: none"> ■ Special anti-fog coating ■ Soft, foam anti-fog strip ■ Soft bindings 	Black/Blue	25/Bx, 4 Bx/Cs
15710	GlareGard™ Procedure Mask with Extended Shield			

10/20/2008

Precept Medical Products, Inc. 370 Airport Road • Arden, NC 28704 • 1-800-438-5827 • www.preceptmed.com