

Testing Documentation

Procedure Masks

Item #	Description	BFE (3 µm)	PFE (0.1 µm)	Delta P	Flammability
15100	Isolation Mask	≥98%	≥98%	Pass	Class 1
15101	Procedure Mask	≥98%	≥98%	Pass	Class 1
15110	Procedure Mask	≥99%	≥98%	Pass	Class 1
15111	Procedure Mask	≥99%	≥98%	Pass	Class 1
15120	Procedure Mask	≥98%	≥97%	Pass	Class 1
15130	KittyKare Procedure Mask	≥99%	≥99%	Pass	Class 1

Surgical Masks

Item #	Description	BFE (3 µm)	PFE (0.1 µm)	Delta P	Flammability
15200	Lite & Cool Surgical Mask	≥98%	≥98%	Pass	Class 1
15201	Standard Surgical Mask - Classic Style	≥98%	≥98%	Pass	Class 1
15215	Sensitive Skin Surgical Mask	≥98%	≥98%	Pass	Class 1
15220	Lite Pouch Surgical Mask	≥99%	≥99%	Pass	Class 1
15230	KittyKare Surgical Mask	≥98%	≥98%	Pass	Class 1
15250	Lite & Cool Horizontal Tie Mask	≥99%	≥98%	Pass	Class 1
15600	Standard Surgical Mask	≥99%	≥98%	Pass	Class 1
65-3333	Comfort-Cone Surgical Mask	≥98%*	N/A	Pass	Class 1

Specialty Masks

*Modified Greene & Vesley Method

Item #	Description	BFE (3 µm)	PFE (0.1 µm)	Delta P	Flammability
15210	Foam Anti-Fog Surgical Mask	≥98%	≥98%	Pass	Class 1
15211	Film Anti-Fog Surgical Mask	≥98%	≥98%	Pass	Class 1
15212	Tape Anti-Fog Surgical Mask	≥98%	≥98%	Pass	Class 1
15216	Sensitive Skin Foam Anti-Fog Surgical Mask	≥98%	≥98%	Pass	Class 1
15525	Laser Mask, Fluid Resistant Anti-Fog	≥99%	≥99%	Pass	Class 1

Fluid Resistant Masks

Item #	Description	ASTM F1862 (Blood Splash Resistance)	BFE (3 µm)	PFE (0.1 µm)	Delta P	Flammability
15300	FluidGard® Procedure Mask	160 mmHg	≥99%	≥99%	Pass	Class 1
15301	FluidGard® Anti-Fog Procedure Mask	160 mmHg	≥99%	≥99%	Pass	Class 1
15310	FluidGard® Procedure Mask w/ Extended Shield	160 mmHg	≥99%	≥99%	Pass	Class 1
15315	FluidGard® Procedure Mask w/ SnapShield	160 mmHg	≥99%	≥99%	Pass	Class 1
15320	FluidGard® Surgical Mask	160 mmHg	≥99%	≥99%	Pass	Class 1
15330	FluidGard® Surgical Mask w/ Extended Shield	160 mmHg	≥99%	≥99%	Pass	Class 1
15335	FluidGard® Surgical Mask w/ SnapShield	160 mmHg	≥99%	≥99%	Pass	Class 1
15340	FluidGard® Procedure VeilMask	120 mmHg	≥99%	≥99%	Pass	Class 1
15345	FluidGard® Procedure VeilMask w/ Extended Shield	120 mmHg	≥99%	≥99%	Pass	Class 1
15350	FluidGard® Surgical VeilMask	120 mmHg	≥99%	≥99%	Pass	Class 1
15355	FluidGard® Surgical VeilMask w/ Extended Shield	120 mmHg	≥99%	≥99%	Pass	Class 1
15400	FluidGard® KittyKare Procedure Mask	160 mmHg	≥99%	≥99%	Pass	Class 1
15410	FluidGard® KittyKare Procedure Mask w/ Extended Shield	160 mmHg	≥99%	≥99%	Pass	Class 1
15420	FluidGard® KittyKare Surgical Mask	160 mmHg	≥99%	≥99%	Pass	Class 1
15500	Fluid Resistant Procedure Mask	160 mmHg	≥99%	≥99%	Pass	Class 1
15510	Fluid Resistant Procedure Mask w/ Extended Shield	160 mmHg	≥99%	≥99%	Pass	Class 1
15520	Fluid Resistant Surgical Mask	160 mmHg	≥99%	≥99%	Pass	Class 1
15530	Fluid Resistant Surgical Mask w/ Extended Shield	160 mmHg	≥99%	≥99%	Pass	Class 1
15710	GlareGard™ Procedure Mask with Extended Shield	160 mmHg	≥99%	≥99%	Pass	Class 1
15730	GlareGard™ Surgical Mask with Extended Shield	160 mmHg	≥99%	≥99%	Pass	Class 1
65-3395	N95 Particulate Respirator (Regular Size)	160 mmHg	≥99%	≥99%	Pass	Class 1
65-3395s	N95 Particulate Respirator (Small Size)	160 mmHg	≥99%	≥99%	Pass	Class 1