

Surgical Face Masks



15201



15220



15250



Precept's surgical masks are designed to offer comprehensive protection while remaining comfortable and easy to don. Bacterial filtration efficiencies for all Precept surgical masks exceed 98%. Each surgical mask incorporates an inner mask material with excellent fiber pindown; this helps prevent facial irritation. Masks are sonically welded, providing for a durable construction. In addition, Precept's surgical masks are cool and breathable, excellent for use during long procedures.

Item #	Description	Features	Color	Packaging
15200	Lite & Cool Surgical Mask	<ul style="list-style-type: none"> Extra soft ties and bindings Keeps wearer cool & comfortable during long procedures 	Blue	50/Bx, 6 Bx/Cs
15201	Standard Surgical Mask (Classic Style)	<ul style="list-style-type: none"> Traditional surgical mask style Excellent fiber pindown helps eliminate facial irritation 	Blue	50/Bx, 6 Bx/Cs
15215	Sensitive Skin Surgical Mask	<ul style="list-style-type: none"> Contains no dyes or colorants that irritate sensitive skin Extra soft ties and bindings Constructed of breathable materials to keep the wearer cool & comfortable 	White	50/Bx, 6 Bx/Cs
15220	Lite Pouch Surgical Mask	<ul style="list-style-type: none"> Pouch style mask keeps material off face to allow for greater comfort and easier breathing 	Blue	50/Bx, 6 Bx/Cs
15230	KittyKare Surgical Mask	<ul style="list-style-type: none"> Extra soft ties and binding Child-friendly face mask with outer print of blue and pink kittens 	KittyKare Print	50/Bx, 6 Bx/Cs
15250	Lite & Cool Horizontal Tie Mask	<ul style="list-style-type: none"> Keeps wearer cool and comfortable during long surgical procedures Excellent bacterial filtration efficiency Horizontal ties 	Blue	50/Bx, 12 Bx/Cs
15600	Standard Surgical Mask	<ul style="list-style-type: none"> Value line mask ideal for surgical and general use procedures Extra soft ties and bindings 	Blue	50/Bx, 6 Bx/Cs

Surgical Face Masks



65-3010



65-3130



65-3342



Precept offers an array of surgical masks, to meet the various needs and challenges of medical personnel. Masks are available in a variety of styles, including traditional, cone, off the face (also called “duckbill” or “vent” style), and with and without extended shields. Precept also offers an exclusive Stretch Knit™ tie that enables the wearer to adjust the tie tension for a snug fit. This helps reduce the “blow-by” that can be experienced when wearing surgical masks. Stretch Knit™ ties stretch and give as the wearer moves his/her head, then return to the previous position, continuing to ensure a proper fit.

Precept’s surgical masks have bacterial filtration efficiencies exceeding 98%, and all masks have Class 1 flammability.

Item #	Description	Features	Color	Packaging
65-3010	Comfort-Cool™ Surgical Mask	<ul style="list-style-type: none"> Constructed of lightweight & comfortable materials Ideal for use during long procedures Enclosed nosepiece for contour fit 	Blue	50/Bx, 6 Bx/Cs
65-3100	Comfort-Soft™ Surgical Mask	<ul style="list-style-type: none"> Constructed of exceptionally soft materials Contour nosepiece provides a snug, comfortable fit 	Blue	50/Bx, 6 Bx/Cs
65-3110	Comfort-Plus™ Surgical Mask	<ul style="list-style-type: none"> Flat fold design to eliminate fabric creases; poly ties Low lint material reduces potential for facial irritation 	Blue	50/Bx, 6 Bx/Cs
65-3112	Comfort-Plus™ Surgical Mask	<ul style="list-style-type: none"> Stretch Knit™ Ties Excellent fiber pindown to reduce facial irritation 	Blue	50/Bx, 5 Bx/Cs
65-3120	Comfort-Plus™ Surgical Mask	<ul style="list-style-type: none"> Flat fold design to eliminate fabric creases; poly ties Low lint material reduces potential for facial irritation 	Green Diamond	50/Bx, 6 Bx/Cs
65-3122	Comfort-Plus™ Surgical Mask	<ul style="list-style-type: none"> Stretch Knit™ Ties Excellent fiber pindown to reduce facial irritation 	Green Diamond	50/Bx, 5 Bx/Cs
65-3130	Comfort-Elite™ Surgical Mask	<ul style="list-style-type: none"> Soflin™ material cushions wear’s sensitive skin at critical contact points Contains no dyes or colorants which irritate skin 	White	50/Bx, 6 Bx/Cs



65-3124



65-3333



65-3348



Item #	Description	Features	Color	Packaging
65-3312	FluidGard™ Mask	<ul style="list-style-type: none"> ■ Stretch Knit™ Ties ■ Enclosed nosepiece 	Blue Diamond	50/Bx, 5 Bx/Cs
65-3314	FluidGard™ Mask	<ul style="list-style-type: none"> ■ Poly ties ■ Enclosed nosepiece 	Blue Diamond	50/Bx, 6 Bx/Cs
65-3333	Comfort-Cone™ Surgical Mask	<ul style="list-style-type: none"> ■ Sturdy nosepiece helps hold mask in place ■ Flanged edge offers comfortable fit 	Blue	50/Bx, 12 Bx/Cs
65-3342	DualGard™ Mask with Extended Shield	<ul style="list-style-type: none"> ■ Snap-on shield provides maximum facial protection ■ Special anti-fog coating applied to shield to ensure visual acuity ■ Stretch Knit™ Ties ■ Black anti-glare strip 	Blue Diamond	25/Bx, 4 Bx/Cs
65-3344	DualGard™ Mask with Extended Shield	<ul style="list-style-type: none"> ■ Snap-on shield provides maximum facial protection ■ Special anti-fog coating applied to shield to ensure visual acuity ■ Cutaway nose section in shield ensures comfortable fit ■ Black anti-glare strip; poly ties 	Blue Diamond	25/Bx, 4 Bx/Cs
65-3346	DualGard™ Vent™ Mask with Extended Shield	<ul style="list-style-type: none"> ■ Vent style mask keeps material off face to allow for greater comfort and easier breathing ■ Cutaway nose section in shield ensures comfortable fit ■ Black anti-glare strip ■ Snap-on shield is fog resistant 	Blue Diamond	25/Bx, 4 Bx/Cs