

## VeilMasks



15350

15355



**Precept** offers an innovative VeilMask design that provides excellent protection against bacteria and splash from bodily fluids. The fluid resistant veil material extends from the bottom of the face mask to the upper chest. This protects the neck area which is left exposed by traditional face masks. The veil can either be worn loose over the gown or tucked into the gown to provide a secure barrier. VeilMasks pass the ASTM F 1862 blood splash test at 120 mm Hg and filter against more than 99% of bacteria as small as 0.1 micron in size. VeilMasks are available with and without extended shields.

Patent No. 5596985

### VeilMask Protection



**Exposed!**



**VeilMask Protection**

The VeilMask offers *ultimate protection* for the exposed neck and chest area!

Item #	Description	Features	Color	Packaging
15340	FluidGard™ Procedure VeilMask	<ul style="list-style-type: none"> <li>■ Sonic welded earloops prevent facial irritation</li> </ul>	Blue Diamond	25/Bx, 10 Bx/Cs
15345	FluidGard™ Procedure VeilMask with Extended Shield	<ul style="list-style-type: none"> <li>■ Extended shield provides maximum facial protection</li> <li>■ Soft foam anti-fog strip prevents fogging of eyeshield</li> <li>■ Black anti-glare strip helps ensure optical clarity</li> </ul>	Blue Diamond	25/Bx, 4 Bx/Cs
15350	FluidGard™ Surgical VeilMask	<ul style="list-style-type: none"> <li>■ Sonic welded construction</li> <li>■ Soft ties and bindings</li> </ul>	Blue Diamond	25/Bx, 6 Bx/Cs
15355	FluidGard™ Surgical VeilMask with Extended Shield	<ul style="list-style-type: none"> <li>■ Extended shield provides maximum facial protection</li> <li>■ Soft foam anti-fog strip prevents fogging of eyeshield</li> <li>■ Black anti-glare strip helps ensure optical clarity</li> </ul>	Blue Diamond	25/Bx, 4 Bx/Cs